DIRECT/NON-COPING TECHNIQUE

Like the original AccessPost, the patented AccessPost Overdenture system is a passive post that provides the necessary strength, retention and stability a restoration requires. It also gives you exactly what the name says: access to the periapical tissues if they become infected after the root canal is completed and the post has long since been cemented. Without sacrificing the post's retention, AccessPost Overdenture gives you the one option you may need to conservatively treat an apical infection - easy post removal without widening the original canal. Retreatment of a failed root canal is possible as a routine and predictable in-office, non-surgical service for every practitioner. This is accomplished by using the special retreatment drills that are available in a separate, optional **AccessPost Retreatment Kit (Cat. No. 570-0-570-04)**. Note: Only the AccessPost retreatment drills should be used. No other drills will perform the task.

1. ENGLISH: Post selection is best done by placing a post in a hemostat and holding it over an undistorted x-ray. If there is a minimum of approximately 0.5mm of lateral root structure at the most apical placement of the post, there is sufficient root structure for placement.

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2. ENGLISH: Use a Peeso or Gates Glidden reamer to remove gutta percha, and for preliminary sizing of the canal (*EDS' Gates Glidden Drills are recommended*). Then use the primary reamer.

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3. ENGLISH: Use the primary reamer to prepare full length of post hole. (*Cut Wet*).

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4. ENGLISH: Use the OVD countersink drill/root facer to create the flange & second tier preparations. (*The flange & second tier must be fully seated. If not, you reduce the post's stability and retention and increase the chance of post fracture under function.*)

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5. ENGLISH: Trial insertion of the post to determine fit.

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
6. ENGLISH: Post Adjustment: Remove post. Cut apical end if necessary to ensure full seating of the flange and second tier. Note: apical adjustment is always done after trial insertion.

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
7. ENGLISH: Post cementation. Note fully seated flange & second tier (A). (*Flexi-Flow Cat. No. 850-00 or Flexi-Flow Natural Cat. No. 860-00 fluoride releasing composite cement from EDS is recommended*). Note: For the highest retention, use the AccessPost Retention Kit Cat.No. 560-00.

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
8. CAP INSERTION
ENGLISH: The dentist is now ready to insert the female cap into the denture. Place nylon female cap on ball of post. Make sure the rubber band is covering the height of contour of the head. If not, there is a risk that the cold cured acrylic could lock in under the head, making removal of the denture difficult.

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
9. ENGLISH: Place marking paste on cap. Seat denture. Relieve denture to make room for cap.

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10. ENGLISH: Place pink cold cure acrylic in relieved denture. Seat denture, let set, and pick up attachment cap.

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
NOTE: The cap should have a little clearance from the coronal tooth structure (or coping) when seated on the ball of the post so that it will rotate. If not, simply relieve a small amount of nylon on lip of cap, and if necessary, relieve some coping or peripheral tooth structure to create space.

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11. ENGLISH: Remove rubber band and trim flash. Overdenture is now ready for use.

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
Refill Instruction Sheet
For Direct Technique and Indirect Technique

AccessPost Overdenture

Caution: Federal law restricts this device to sale by or on the order of a dentist.
NON-STERILE

STAINLESS STEEL POSTS
PROFFEN AUS ROSTFREIEM STAHL
PIVOTS EN ACIER INOXYDABLE
PERNI IN ACCIAIO INOSSIDABILE
ESPIGAS DE ACERO INOXIDABLE

INDIRECT/COPING TECHNIQUE

1 - 3. ENGLISH: Follow steps 1 - 3 of "direct/non-coping technique" (see reverse side).

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4. ENGLISH: Use the secondary drill to prepare seat for second tier.

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
5. ENGLISH: Trial insert the post to determine fit, making sure that the second tier is fully seated. The flange of post must sit on the coronal tooth structure.

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6. ENGLISH: Post adjustment: Remove post. Cut apical end if necessary to ensure seating of second tier and coronal contact of flange. Note: apical adjustment is always done after trial insertion.

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
7. 0.25mm ENGLISH: Seat the corresponding size brass transfer stud (#of rings = size of post) allowing for a 0.25mm space between the flange and the coronal tooth structure for the coping to be placed.

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
8. ENGLISH: Take an impression of the seated brass transfer stud.

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9. ENGLISH: Remove the impression (including the brass transfer stud from the canal) and temporarily seal the canal.

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ENGLISH: Remove the impression (including the brass transfer stud from the canal) and temporarily seal the canal.

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10. ENGLISH: Send the impression to the laboratory along with the appropriate size brass transfer stud. (**#of rings = size of post**).

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
11. ENGLISH: The laboratory lubricates and pours it up in stone or plaster.

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

12. ENGLISH: The laboratory then waxes up and casts the coping.

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13. ENGLISH: Note: When cementing AccessPost Overdenture, A. Use composite cement only in a canal that has been either etched or grooved (for canal grooving, EDS recommends EDS grooving drills Cat. No. 593-0 & Cat. No. 593-123) and has had all remnants of any eugenol containing temporary cements or root canal pastes removed. B. Do NOT use glass ionomer cements, due to their brittle nature. The dentist cements the coping (a). While the cement is wet AccessPost Overdenture post is cemented in place (b). (Flexi-Flow Cat. No. 850-00 or Flexi-Flow Natural Cat. No. 860-00 fluoride releasing composite cement from EDS is recommended). **Note: For the highest retention, use the AccessPost Retention Kit Cat. No. 560-00.** The insertion stops when the flange is fully seated within the coping.

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14. ENGLISH: Follow instructions 8 - 11 of "direct/non-coping technique" for cap incorporation.

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For the latest product information, please visit www.edsdental.com

89 Leaning Street, S. Hackensack, NJ 07606, U.S.A.

201-487-9090 • info@edsdental.com

MADE IN USA

U.S. Patent No. 5,275,563 & 5,413,480 and Foreign Patents Pending.