IMPORTANT:

Read pages 7 through 12 for
Technique first.

Caution: Federal law restricts this device to sale by or on the order of a licensed dentist.

Flexi-Flange[®] Fiber


*Instruction Book for Flexi-Flange
Fiber Posts*

Recommended uses for Flexi-Flange Fiber

| Post Number | 0 | 1 | 2 |
|-------------|--|---|---|
| Color Code | Yellow | Red | Blue |
| | thin to average buccal or mesial roots of molars | average to large buccal or mesial roots of molars | average roots of all maxillary anteriors |
| | thin to average roots of maxillary first premolars | normal to large roots of maxillary premolars | average to large roots of premolars |
| | average roots of mandibular incisors | average roots of maxillary laterals | large roots of mandibular anteriors |
| | | average distal and palatal roots of molars | average to large distal and palatal roots of molars |
| | | | average maxillary centrals |
| | | | average maxillary canines |

Note: In all Flexi-Flange Fiber restorations, a minimum of a 1.5mm ferrule around the tooth is needed for the final restoration.


Technique: Use of the Depth Gauge in Post Selection

Research shows that parallel, solid shanked posts should have at least 1 mm of tooth structure lateral to the most apical placement of the post. To aid in this placement, Flexi-Flange Fiber uses a transparent plastic depth gauge with silhouettes of the different sizes of the posts. Lateral and parallel to each silhouette are vertical lines spaced 1 mm from the threads. By placing the gauge over an accurate radiograph of a tooth, the dentist may easily determine if the 1 mm of lateral clearance exists. If the lines fall outside the root on the x-ray there is potentially not enough lateral tooth structure for safe placement. In the latter case, the dentist should either go to a smaller post or remove enough apical post length for the post to fit at least 1 mm within the external borders of the root. Most importantly, the flange of the post must always be fully seated. The dentist should never allow the flange of the post not to seat. This would allow a loose coronal fit that would increase the chances of the post loosening over time!

If the dentist chooses to remove apical length of the post, either because the full length of the placed post would thin out the lateral tooth structure too much or because the post-hole is too short for placement of the complete post length, he should follow the steps listed below:

- 1) Identify which post size to place by using the depth gauge placed over an xray.
- 2) Choose the correct tap size that corresponds to the post size
- 3) Thread the tap into the post-hole thus creating the internal thread in the root for the post. Do not force the tap. If it does not fully seat back up a 1/4 turn and then proceed forward.
- 4) Unthread the tap from the root.
- 5) Trial seat the post
- 6) Cut off the necessary apical post length (Do not cut tap), allowing the flange to seat fully.
- 7) Cement the post.



Post-Hole Preparation

The post-hole preparation begins with the removal of the root filling material using either a Peeso or Gates Glidden reamer. Then, in sequence, a non-end cutting drill (Peeso or Gates Glidden reamer) is used until 100% of the post-hole length and 90% of the post-hole width have been established. The following chart indicates which non-end cutting drill will produce 90% of the post-hole width for the various Flexi-Flange Fiber sizes.

| Peeso | | Gates Glidden | | EDS Gates Glidden | | Flexi-Flange Fiber Primary Reamer |
|-------|----|---------------|----|-------------------|--------|-----------------------------------|
| 1 | or | 2 | or | yellow | then ↘ | 0 (yellow) |
| 3 | or | 4 | or | red | then ↘ | 1 (red) |
| 4 | or | 5 | or | blue | then ↘ | 2 (blue) |

When 100% of the post-hole length and 90% of the width have been achieved, the primary reamer is used. Since the Flexi-Flange Fiber will fit optimally if a more concentric hole is maintained, **the number of entries into the post-hole with the primary reamer should be limited**. It is much easier to prepare the post-hole when the canal is lubricated with either water or an anesthetic solution, or with any suitable wetting agent.

The countersink drill cuts two preparation in one operation. It prepares the seat for the second tier, as well as the seat for the flange of the post. The post must always fully seat. You can determine full seating of the post by the flush fit of the flange within its preparation. If the coronal flat surface of the root is slanted to the buccal, the flange may be seated lingually, but not buccally. In this case, deepen the countersink preparation enough to assure the flush seating of the buccal portion of the post. There is no danger in countersinking the post too much. If on the other hand, the dentist does not seat the post completely, he is reducing Flexi-Flange Fiber's retention and increasing the chances of the post loosening or fracturing under function. To achieve a complete seating, in post-hole preparations shorter than the length of the shank of the post to be placed, the dentist must remove enough apical post length to allow full seating of the post's second tier and flange. The countersink drill is now used to prepare the post-hole for the second tier and flange portions of the Flexi-Flange Fiber. The smooth extension on the drill is simply a lead to facilitate parallelism between the primary post-hole, second tier and flange.


Tap Insertion

Note: For optimum retention it is recommended that you flush the canal with 17% EDTA for 1 minute before and during tapping the canal. (EDS recommends the use of its EDTA Plus, Cat. No. 770-16). This will remove debris and open the dentinal tubules, thereby increasing post retention.



After choosing the correct size post with the depth gauge, the correct size tap that corresponds with that post is inserted. It is important to note that the tap is designed to be seated with light pressure. During the tap-seating, if moderate resistance is felt, the tap may be backed off 1/4 to 1/2 turn and then advanced again. Advancing while backing off 1/4 turn when moderate resistance occurs is repeated until the tap is fully inserted and the thread is created inside the root canal for the post. This procedure will remove debris from the thread line and facilitate insertion. NOTE: Please use appropriate autoclave sterilization methods such as heat or steam prior to each use. Additionally, periodic ultrasonic cleaning to remove debris from the split-shank is recommended.

***Important Note:** If the primary reamer did not reach its' full depth the tap will not seat fully. Tap the root to the length of the primary reamer and shorten the post apically so that the flange of the post is fully seated.

Flush the canal with 17% EDTA

for 1 min.


Fully Seat Tap
(See Important Note)

Tap Insertion while wet with 17% EDTA

Post Insertion

The tap creates the thread pattern inside the root canal for the post. The correct size post is now ready to be trial seated. Using the Flexi-Flange Fiber wrench thread post with light pressure until it is fully seated.

The post is now unthreaded out of the root. At this point, alterations to the post may be made. It is extremely important to note that the flange must always fully seat. Therefore, alteration should be made to the apical end of the post. If needed, cut the apical end with a diamond disk until the flange is fully seated.


*Post Adjustment
(If necessary)*

Cementation

Dry the canal


Place composite resin cement in the canal - **No Bonding Agent is Required**. Place post on wrench and thread until fully seated. Remove the wrench and follow with core build-up and preparation. Please note: For the greatest post retention we recommend the use of Flexi-Flow Cem[®] (Cat. No. 850-00) or, Flexi-Flow Natural[®] (Cat. No. 860-00) or, Flexi-Flow Auto[®] (Cat. No. 870-00) or, Flexi-Flow Auto[®] E (Cat. No. 880-00) fluoridated composite cements. **Due to the threads of the Flexi-Flange Fiber post, post retention is greater than passive posts cemented with a bonding agent. Therefore, a bonding agent is not needed.**

The post is inserted into the post-hole and threaded in with light pressure. The post will seat completely with minimal resistance. Excess cement is now removed. The Flexi-Flange Fiber has now been inserted and cemented with minimal stress being transmitted to the root.

Core Formation

Composite Cores - EDS recommends the use of clear core forms (Cat. No. 890-00) to be used with composite materials. The composite should be placed in the core form and seated over the post, using moderate pressure to ensure close adaptation of the composite to the core. For best results Essential Dental Systems recommends the use of Ti-Core® (Chemical Cured, Cat. No. 800-00) or, Ti-Core® Natural (Chemical Cured, Cat. No. 810-00), Ti-Core® Auto E (Dual Cured, Cat. No. 830-00) or, Ti-Core® Auto White (Dual Cured, Cat. No. 840-00) fluoridated composite core materials.

For more information and the latest instructions please visit
www.edsdental.com

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1-800-22-FLEXI

Or contact your

Essential Dental Systems, Inc. dealer.

Essential Dental Systems, Inc.

89 Leuning Street, S. Hackensack, NJ 07606


Tel: 201-487-9090 Fax: 201-487-5120

www.edsdental.com

Flexi-Flange is a trademark of Essential Dental Systems, Inc.

© Essential Dental Systems, Inc.

Made in U.S.A.

WARNING: This product can expose you to the chemical Bisphenol A, which is known to the State of California to cause reproductive harm. For more information go to www.P65Warnings.ca.gov.

ENDODONTIC FIBER POST CLEANING

While chemical sterilants such as peracetic acid, 2% glutaraldehyde based formulations, and/or 6% hydrogen peroxide may be used for cleaning endodontic fiber posts, the majority of fiber endodontic post manufacturers recommend the use of disinfectants such as 70% ethyl or isopropyl alcohol^{1,2,3,4,8}. This has become a dental industry standard for the cleaning of all fiber post systems before placement. The following protocol is recommended by most dental industry's top manufacturers including 3M, Brassler, as well as Dentsply Sirona^{5,6,7}:

- Disinfect the post with alcohol and dry it with air free of water and oil.
- Wipe the post with alcohol.
- Allow the post to air dry for 30 seconds before use.

While not necessary, Endodontic glass fiber posts may also be sterilized using conventional methods according to a study by *Yagci et al*⁴ without any deleterious structural effects to the post. In this study, conventional steam autoclave sterilization was performed on a large number of glass fiber posts and then the bond strength, flexural strength as well as the elastic modulus of each was evaluated. Additionally, the microstructures of the tested posts were evaluated by SEM and energy-dispersive X-ray spectroscopy (EDS) for damage resulting from the autoclave steam sterilization process. *One-way ANOVA indicated no statistically significant differences* among the groups means in terms of bond strength ($p > 0.05$), flexural strength ($p > 0.05$), or elastic modulus ($p > 0.05$). The results of this study show that glass-fiber posts can be safely and effectively sterilized either by autoclave when necessary, without any negative effect on bond strength, flexural strength, or elastic modulus⁴.

While studies have shown that glass fiber posts are not negatively affected during the steam autoclave sterilization process^{2,3,4,8}, cleaning of glass fiber posts with alcohol ($\geq 70\%$) is adequate for proper disinfection and recommended as per industry standard protocols^{5,6,7}. The cleaning, disinfection, and sterilization of glass fiber posts is ultimately left to the discretion of the clinician as dental industry standard protocol and/or steam autoclave sterilization procedures are both effective/safe. The disinfection of glass fiber posts using alcohol has long been an industry standard with an extensive history of safe clinical use and efficacy.

Citations

1. Bateman G, Ricketts DNJ, Saunders WP. Fiber-based post systems: a review. Br Dent J 2003;195:43–48.
2. Farrugia C, Cassar G, Valdramidis V, Camilleri J. Effect of sterilization techniques prior to antimicrobial testing on physical properties of dental restorative materials. J Dent 2015;43:703–714.
3. Kotha SB, Ramakrishnaiah R, Devang DD, Celur SL, Qasim S, Matinlinna JP. Effect of disinfection and sterilization on the tensile strength, surface roughness, and wettability of elastomers. J Investig Clin Dent 2016;8: doi:10.1111/jicd.12244.
4. Yagci F, Ustun Y, Zortuk M, Agirnasligil M. Effect of Sterilization on Bond Strength and Mechanical Properties of Fiber Posts. J Adhes Dent. 2019;21(2):143-148. doi: 10.3290/j.jad.a42325. PMID: 30949627.
5. <https://multimedia.3m.com/mws/media/3858070/3m-relyx-fiber-post-cementing-endodontic-post.pdf>.
6. <http://brasslerusa.com/wp-content/files/IFU-0026%20ES%20Fiber%20Post%20System%20IFU%20REV%20A.pdf>
7. <https://www.dentsplysirona.com/content/dam/master/product-procedure-brand-categories/endodontics/product-categories/ifu/END-IFU-Endodontic-Sterilization-Procedure-EN-1402.pdf>
8. Patel SS, Sethuraman R. A randomized controlled twelve month clinical study on the evaluation of success rate of endodontically treated teeth restored with metal poly-fiber posts and dentin posts. J Indian Prosthodont Soc. 2022 Jan-Mar;22(1):38-47. doi: 10.4103/jips.jips_134_21. PMID: 36510946; PMCID: PMC8884352.