## Ti-Core Auto E Instrucciones

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## Ti-Core Auto E Instructions

Ti-Core® Auto E is a dual cured hybrid composite core material reinforced with the lanthanide series of elements.


Ti-Core Auto E is formulated to match Shade Vita A2.

### AREA OF USE:

- 1. NONVITAL TEETH** - Ti-Core Auto E is recommended for use with all prefabricated posts and/or pin build-ups.
- 2. VITAL TEETH** - Ti-Core Auto E may be used to build-up missing tooth structure. A calcium hydroxide or glass ionomer base may be used under the Ti-Core Auto E composite.
- 3. BONDING TO ENAMEL & DENTIN:**
  - a. Ti-Core Auto E can be bonded to enamel with the acid etch technique.
  - b. Ti-Core Auto E can be bonded to dentin using any of the available multi-step dentin bonding agents currently on the market.

### Clinical Technique:

1. Remove the syringe cap, place the mixing stator on the syringe and place the delivery tip onto the mixing stator. (Figures 1-4)

- 1.** Before attaching the mixing tip, remove the syringe cap by turning 90° counter clockwise and pulling to remove.
- 2.** Push the mixing tip onto the syringe. Line up the notch on the tip with the corresponding notch on the syringe.
- 3.** Turn the tip 90° clockwise until it locks into position.
- 4.** Remove and discard the mixing tip after use and replace the syringe cap.

2. As pressure is applied to the plunger, the base and catalyst are extruded through the stator and mixed automatically. This guarantees a consistent mix each and every time.

### Procedural Timing:

Working Time: 3:45 minutes

Light Cure: 20 seconds (< 2mm)

40 seconds (< 4mm)

Chemical Cure: 6 minutes intra oral

Ti-Core Auto E has a 3:45 minute working time. The material may be manipulated using a composite instrument. A bonding agent can be used as a separating medium to allow for easier manipulation.

**Chemical Cure:** When using only the auto cure mode Ti-Core Auto E has a 6 minute setting time (37°C/98°F). Lower temperatures delay the set.

**Light Cure:** The material may be light cured with the usual dental halogen lamps for at least 40 seconds for 4mm maximum layer, anytime after extrusion from the syringe. A thin layer of less than 2mm will need only a 20 second exposure.

### Physical Properties

Working Time: 3:45 minutes

Setting Time: 6:00 minutes

% Filled Approx. 58%

Range of partial size Sub-micron to 10 micron

Store between 60°F (15°C) - 85°F (29°C).

Store unopened boxes in the refrigerator for maximum shelf life.

Do not freeze.

Use at room temperature to ensure proper set time!

Do not expose to direct sunlight.

Do not use after expiration date.

MSDS sheets available. Call: 1-800-22-FLEXI or 201-487-9090. E-mail: info@edsdental.com.

### Precautions:

- If close to pulp use a suitable liner.
- Contact with eyes and skin should be avoided. If accidental contact occurs rinse immediately with water and consult a physician if necessary.
- Do not use materials containing eugenol it will delay or inhibit the set.
- Water and oily air inhibit polymerization.
- A dry field is indicated when placing a composite.
- There are no known systemic side effects to date, however, do not use in patients with known sensitivities to Bis-GMA, benzyl peroxide or composite resin materials.


# Ti-Core® Auto E

## Dual Cured Reinforced Core Material

Cat. No. 830-00

Tooth Color Vita® A2

Caution: Federal law restricts this device to sale by or on the order of a licensed dentist.

Releasing Fluoride 10+ Years & Counting

**EDS** ESSENTIAL DENTAL SYSTEMS®

89 Leuning Street  
S. Hackensack, NJ 07606  
Phone: 1-800-22-FLEXI • 201-487-9090

Fax: 201-487-5120

www.edsdental.com

U.S. Patent Nos. 5,034,433 and 5,204,398.  
European Patent No. EP 0 631 598 B1. Other patents pending.

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