## Ti-Core Auto E Instrucciones

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## Ti-Core Auto E Instrucciones

Ti-Core® Auto E is a dual cured hybrid composite core material reinforced with the lanthanide series of elements.


Ti-Core Auto E is formulated to match Shade Vita A2.

### AREA OF USE:

- 1. NONVITAL TEETH** - Ti-Core Auto E is recommended for use with all prefabricated posts and/or pin build-ups.
- 2. VITAL TEETH** - Ti-Core Auto E may be used to build-up missing tooth structure. A calcium hydroxide or glass ionomer base may be used under the Ti-Core Auto E composite.
- 3. BONDING TO ENAMEL & DENTIN:**
  - a. Ti-Core Auto E can be bonded to enamel with the acid etch technique.
  - b. Ti-Core Auto E can be bonded to dentin using any of the available multi-step dentin bonding agents currently on the market.

### Clinical Technique:

1. Remove the syringe cap, place the mixing stator on the syringe and place the delivery tip onto the mixing stator. (Figures 1-4)

- 1.** Before attaching the mixing tip, remove the syringe cap by turning 90° counter clockwise and pulling to remove.
- 2.** Push the mixing tip onto the syringe. Line up the notch on the tip with the corresponding notch on the syringe.
- 3.** Turn the tip 90° clockwise until it locks into position.
- 4.** Remove and discard the mixing tip after use and replace the syringe cap.

2. As pressure is applied to the plunger, the base and catalyst are extruded through the stator and mixed automatically. This guarantees a consistent mix each and every time.

Procedural Timing:		
Working Time:	Approx. 3 minutes	Ti-Core Auto E has an approximate 3 minute working time. The material may be manipulated using a composite instrument. A bonding agent can be used as a separating medium to allow for easier manipulation.
Light Cure:	20 seconds (< 2mm)	40 seconds (< 4mm)
Chemical Cure:	Approx. 5 minutes intra oral	

**Chemical Cure:** When using only the auto cure mode Ti-Core Auto E has an approximate 5 minute setting time (37°C/98°F). Lower temperatures delay the set.

**Light Cure:** The material may be light cured with the usual dental halogen lamps for at least 40 seconds for 4mm maximum layer, anytime after extrusion from the syringe. A thin layer of less than 2mm will need only a 20 second exposure.

**Physical Properties**

Working Time:	Approx. 3 minutes
Setting Time:	Approx. 5 minutes
% Filled	Approx. 58%
Range of partial size	Sub-micron to 10 micron
Store between 60°F (15°C) - 85°F (29°C).	
Store unopened boxes in the refrigerator for maximum shelf life.	
Do not freeze.	
Use at room temperature to ensure proper set time!	
Do not expose to direct sunlight.	
Do not use after expiration date.	
MSDS/SDS sheets available. Call: 1-800-22-FLEXI or 201-487-9090. E-mail: info@edsdental.com.	

- Precautions:**
- If close to pulp use a suitable liner.
  - Contact with eyes and skin should be avoided. If accidental contact occurs rinse immediately with water and consult a physician if necessary.
  - Do not use materials containing eugenol it will delay or inhibit the set.
  - Water and oily air inhibit polymerization.
  - A dry field is indicated when placing a composite.
  - There are no known systemic side effects to date, however, do not use in patients with known sensitivities to Bis-GMA, benzyl peroxide or composite resin materials.

Please Note:  
Mixing tips and mixing pad sheets are single use only. If necessary, the exterior of the syringe may be disinfected with alcohol or sodium hypochlorite. Do NOT autoclave.

To download Material Safety Data Sheet (MSDS) / SDS (Safety Data Sheets)  
Please Visit: [Http://edsdental.com/sds](http://edsdental.com/sds)

# Ti-Core® Auto E

## Dual Cured Reinforced Core Material

REF 830-00 Tooth Color Vita® A2

Contains: 2 auto-mix syringes (containing 9gm each), 20 mixing tips, 20 intra oral tips.

# Ti-Core® Auto E

## Dual Cured Reinforced Core Material

REF 830-00 Tooth Color Vita® A2

**STORE BETWEEN**  
60°-85°F /  
15°-29°C

**EC REP**  
Authorized European Representative  
European Healthcare & Devices Ltd.  
Stratton House  
Bishopstown Road  
Cork  
T12 Y9TC, Ireland  
Email: [TQMUK@aol.com](mailto:TQMUK@aol.com)

**CE**  
2797

**EDS** ESSENTIAL DENTAL SYSTEMS  
89 Leuning Street  
S. Hackensack, NJ 07606  
Phone: 1-800-22-FLEXI • 201-487-9090  
Fax: 201-487-5120  
[www.edsdental.com](http://www.edsdental.com)

Caution: Federal law restricts this device to sale by or on the order of a licensed dentist.

60 F 15 C 85 F 29 C Caution, Consult Instructions for Use

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