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**ENGLISH: Ti-Core Flow+ Instructions**

Indications for use: Ti-Core Flow+ may be used as a Core Material (use wider syringe tip) or Post Cement (use narrow/longer syringe tip). Please see appropriate section for directions. **Contra-indications:** Do not use in patients with known sensitivities to benzyl peroxide or methacrylate resin materials.

**Procedural Timing:**

Working Time: Approx. 3 minutes Light Cure: 20 seconds (< 2mm) Chemical Cure: Approx. 5 minutes intra oral 40 seconds (< 4mm)

**Chemical Cure:** When using only the auto cure mode Ti-Core Flow+ has an approximate 5 minute setting time (37°C/98°F). Lower temperatures delay the set.

**Light Cure:** The material may be light cured for at least 40 seconds for 4mm maximum layer, anytime after extrusion from the syringe. A thin layer of less than 2mm will need only a 20 second exposure.

**Physical Properties**


% Filled Aprox. 58%  
Range of particle size Sub-micron to 10 micron  
Store between 60°F (15°C) - 85°F (29°C).  
Store unopened boxes in the refrigerator for maximum shelf life.  
Do not freeze.  
Use at room temperature to ensure proper set time!  
Do not expose to direct sunlight.  
Do not use after expiration date.  
MSDS/SDS sheets available. Call: 1-800-22-FLEXI or 201-487-9090. E-mail: info@edsdental.com.

**Precautions:**



- Certain light cured bonding agents may accelerate curing. This may lead to a shorter working time during post cementation.
- If close to pulp use a suitable liner.
- Contact with eyes and skin should be avoided. If accidental contact occurs rinse immediately with water and consult a physician if necessary.
- Do not use materials containing eugenol, it will delay or inhibit the set.
- Water and oily air inhibit polymerization.
- A dry field is indicated when placing a composite.
- Due to the amount of radiopacity, it is recommended to take other precautions to detect voids and secondary caries.

**AREA OF USE:**

- NONVITAL TEETH** - Ti-Core Flow+ is recommended for use with all prefabricated posts and/or pin build-ups.
- VITAL TEETH** - Ti-Core Flow+ may be used to build-up missing tooth structure. A calcium hydroxide or glass ionomer base may be used under the Ti-Core Flow+ composite.
- BONDING TO ENAMEL & DENTIN:**
  - Ti-Core Flow+ can be bonded to enamel with the acid etch technique.
  - Ti-Core Flow+ can be bonded to dentin using any of the currently available multi-step dentin bonding agents on the market.

**Core Material Instructions****Clinical Technique:**

- Remove the syringe cap, place the mixing stator on the syringe and place the delivery tip onto the mixing stator.
- As pressure is applied to the plunger, the base and catalyst are extruded through the stator and mixed automatically. This guarantees a consistent mix each and every time.
- Remove and discard the mixing tip after use and replace the syringe cap.

**Fig 1****Fig 2****Post Cement Instructions****Clinical Technique:**

- Prepare the post-hole in order to receive the cement.
- Remove syringe cap, place the mixing stator on the syringe and place the post-hole delivery tip onto the mixing stator.
- As pressure is applied to the plunger the base and catalyst are extruded through the stator and mixed automatically. This guarantees a consistent mix each and every time.
- Place the post-hole delivery tip into the post-hole and express Ti-Core Flow+ directly into the post-hole (Figure #1).
- In order to achieve maximum penetration depth of the cement into the preparation, use a lentulo spiral to spin the cement to the apical end of the post-hole (Figure #2).
- The post is now ready for core build-up.

Note: Remove and discard the mixing tip after use and replace the syringe cap.  
For maximum tensile strength, the post-hole should be etched.

**Cementation Instructions For Crowns, Bridges, Inlays and Onlays**


- Clean the preparation/cavity with water and then dry excess water using a gentle stream of air for 2 seconds. Do not over dry the dentin.
- Etch the preparation/cavity with 35% phosphoric acid, according to the manufacturer's instructions.
- Condition the contact surfaces of the restoration according to the manufacturer's instructions for use.
- Remove the syringe cap, place the mixing stator on the syringe and place the delivery tip onto the mixing stator.
- As pressure is applied to the plunger, the base and catalyst are extruded through the stator and mixed automatically. This guarantees a consistent mix each and every time.
- Express Ti-Core Flow+ directly into the inner surfaces of the restoration and/or if necessary (to prevent trapped air voids with concave shapes) to the preparation.
- Place the restoration in position using slight pressure.
- Remove rough excess with a brush or spatula and hold the restoration in position with increased pressure. **Note: Remove and discard the mixing tip after use and replace the syringe cap.**

**Ti-Core® Flow+**Combined Dual Cured Reinforced Core Material  
& Endodontic Post CementREF 832-00  
Cat. No.**EDS** ESSENTIAL  
DENTAL  
SYSTEMS

Contains: 2 auto-mix syringes (each containing 4.5 gm base, 4.5 gm catalyst), 20 mixing tips, 20 intra oral tips.  
U.S. Patent Nos. 5,034,433 and 5,204,398. European Patent No. EP 0 631 598 B1. Other patents pending.

**Ti-Core® Flow+**  
Combined Dual Cured Reinforced Core Material  
& Endodontic Post Cement**Shade Vita A2®****EDS** ESSENTIAL  
DENTAL  
SYSTEMS89 Leuning Street  
S. Hackensack, NJ 07606  
www.edsdental.com

Caution: Federal law restricts this device to sale by or  
on the order of a licensed dentist.  
U.S. Patent Nos. 5,034,433 and 5,204,398.  
European Patent No. EP 0 631 598 B1.  
Other patents pending.

TCPSL-4

**English****Please Note:**

Mixing tips and mixing pad sheets are single use only. If necessary, the exterior of the syringe may be disinfected with alcohol or sodium hypochlorite. Do NOT autoclave.  
To download Material Safety Data Sheet (MSDS) / SDS (Safety Data Sheets). Please Visit: [Http://edsdental.com/sds](http://edsdental.com/sds)

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