



EZ-Fill Xpress

Endodontic Obturation System



Notas del Fabricante:

1. La espiral bi-direccional está disponible en acero y níquel-titanio.
2. Es esencial usar el cemento de resina epoxy EZ-Fill diseñado especialmente para la espiral bi-direccional. La resina epoxy EZ-Fill ha sido especialmente formulada para ser radiopaca y tener una viscosidad ideal.
3. El tope de goma azul identifica la espiral bi-direccional de acero. El tope de goma rojo identifica la espiral bi-direccional de Níquel-Titanio.

Para más información y últimas instrucciones, por favor, visiten www.edsdental.com

For optimum obturation Essential Dental Systems recommends the use of the Endo-Express® / SafeSiders® Instrumentation system. (EDS Intro Kit Cat. No. 5600-21, 5600-25, 5610-21, 5610-25).

EZ-Fill Xpress Introductory Kits:

Stainless Steel (Cat. No. 1620-00) & Nickel Titanium (Cat. No. 1625-00).
Each introductory kit contains:

4 - Color coded to ISO Size 25 EZ-Fill bi-directional spirals (1- 21mm length & 3- 25mm length), Two dual barrel syringes (9.5gm each), 20 mixing tips.

EZ-Fill Xpress Refill Kit:

1628-00 EZ-Fill Xpress Epoxy Root Canal Cement
Two - Dual barrel syringe (9.5gm each), 20 mixing tips.

EZ-Fill Bi-Directional Spiral Refill Kits:

Stainless Steel:
1600-01 4-Bi-Directional Spirals (1-21mm length & 3-25mm length)
1600-21 4-Bi-Directional Spirals (21mm length)
1600-25 4-Bi-Directional Spirals (25mm length)

Nickel Titanium:
1605-01 4-Bi-Directional Spirals (1-21mm length & 3-25mm length)
1605-21 4-Bi-Directional Spirals (21mm length)
1605-25 4-Bi-Directional Spirals (25mm length)

EC REP

Authorized European Representative:
Regulatory Affairs ONLY:
European Device Solutions
33, Stapleford Road,
Middlesbrough, Cleveland,
United Kingdom.
Contact E-Mail: TQMUK@aol.com



© 8/06 DL, EZSI-1

89 Leuning Street
S. Hackensack, NJ 07606
Phone: 1-800-22-FLEXI • 201-487-9090
Fax: 201-487-5120
www.edsdental.com