

## Endo-Express™ Mini Latch Head Contra Angle (E-Type)

Standard  
Contra Angle

Mini Head  
Contra Angle

### FEATURES

10 mm Head

E-Type Attachment

Miniature Head Contra  
Angle

### BENEFITS

Allows 22 mm of the gates glidden length to exit the head opening versus 18 mm for a standard contra angle.

Allows for teeth up to 25 mm in length to be prepared to within 3 mm of the apex with the No. 2 gates glidden.

Allows 23 mm of length of the Peeso to exit the contra angle opening allowing it to be used in canals as long as 29 mm length.

In mouths with little access, the mini head often allows the use of shorter shanked instruments. This allows for easier vertical placement of the shorter length.

It eliminates the need to have two slow speed contra angles

It fits onto the same slow speed E-Type motor as your Endo-Express handpiece.

E-Type 1:1, 30,000 RPM, Miniature Latch Head

### Part of the SafeSiders & Endo-Express systems:

These shaping systems allow canals to be shaped in 2-5 minutes.

Replace instruments because they dull, not because you fear they will break!

Reciprocation creates a balanced force, the instruments stay centered.

The 45 degree motion virtually eliminates:

- cyclic fatigue
- torsional stress
- binding
- which virtually eliminates instrument separation

Act as chisels in the clockwise and counter clockwise motion allowing you to go through the tooth and remove debris with unexpected ease.

Prevents canal distortion.

Costs 90% less then rotary NiTi on a per use basis.