



# Essential Organizers

Unique storage systems that let you keep your Essential Dental System posts and accessories organized and sterilized.

## FEATURES / BENEFITS

1. Perfect organizer for office or individual operatories. Keeps posts neat and accessible.
2. Cost effective introductory kit in full or mini size.
3. Convenient infection control: autoclavable, chemclavable or dry heat sterilization.
4. Durable aluminum construction for long-term use.
5. Linear organization of posts for easy identification of correct size and related accessories.



## TECHNICAL FACTS:

Dimensions:

Large Sterilization Box -  
4 1/2" (length), 1 3/4" (width),  
1 1/2" (depth)

Mini-Steri Box - 2" (length),  
1 5/8" (width), 1 1/2" (depth)

## ORDERING INFORMATION:

	Catalog Number	
Sterilization Box (Empty) For storage of 50 posts - 10 posts ea. of 5 sizes (00-03) and accessories	197-00	
Mini-Sterilization Box (Empty) For storage of 10 posts - 2 posts ea. of sizes (00-03) and accessories	197-04	
Flexi-Post Sterilization Box (Complete) Contains 50 posts - 10 posts ea. of sizes (00-03) & accessories	<u>Stainless Steel</u> 197-01	<u>Titanium</u> 197-02
Flexi-Post Mini-Sterilization Box (Complete) Contains 10 posts - 2 posts ea. of 5 sizes (00-03) & accessories	197-05	197-06
Combi-Sterilization Box Contains 16 posts - 12 Flexi-Posts (4 posts ea. of 3 Sizes/0,1,2) and 4 Overdenture posts (2 posts ea. of 2 sizes/1,2) and accessories	197-07	197-03
Flexi-Flange Mini-Sterilization Box (Complete) Contains 8 posts - 2 posts ea. of 4 sizes (0-03) & accessories	197-08	197-09
Special order Combi-Boxes Call 1-800-22-FLEXI to special order.		