



SafeSiders® HF Helical Flat Rotary Instruments

Endo-CSV™ (Control, Shape, Vision) Cordless Endo Motor



**The Affordable
Endo-CSV Cordless
Endodontic Motor
with Built-In LED
Illumination!**

Feature

**Innovative
Endodontist
Designed and
Tested:**

Benefits

Patent pending **relieved helical flat** revolves in the opposite direction of the cutting flutes:

- Increases the flexibility of the instrument
- Reduces engagement along length
- **Creates easier negotiation of canals**
- Produces faster instrumentation to the apex
- Creates less resistance
- Provides a coronal pathway for debris
- Minimizes apical extrusion of debris

Proprietary, NiTi instrument manufacturing process coupled with the relieved helical flat produces:

- An instrument **highly resistant to breakage**
- Instrumentation flexibility which follows the anatomy of the canal, that reduces the risk of ledging, transportation, perforation and minimizing the loss of dentin.
- A system with a reduced number of instruments required for cleaning and shaping while preserving the anatomy.

Economical:

- Time saving rotary instrument shaping system.
- Tough and reusable. SafeSiders HF instruments can be used up to 3 times, providing dramatic savings compared to single use systems*.
- Affordable SafeSiders HF instruments can be used with your current endo motor or, for optimal use and safety the EDS Endo-CSV Cordless Endodontic Handpiece

* Instruments are recommended for 3 time use only in the Endo-CSV Endo Motor. One time use is recommended for extreme curves and use with other endo motors.

SafeSiders HF & Endo-CSV Motor

- Cordless Endodontic Handpiece with Built-in LED Illumination.
- Convenient and completely portable.
- Easy to read LCD display.
- Heavy-duty high torque motor.
- Automatic forward/reverse functions.
- Switch from rotary to reciprocating with a press of a button.
- Five torque settings and nine rotation speeds.
- High capacity lithium-ion battery.
- Operator friendly, easy programming



ORDERING INFORMATION

Endo-CSV Cordless Endodontic Handpiece:

- | | |
|---------|--|
| 5121-00 | Endo-CSV Handpiece, SafeSider HF Instruments in 21 mm length (2 Glidepath Instruments, 2 Shapers, 2 Finishers) 1 box 25/04 GP points, 1 box 30/04 GP points |
| 5125-00 | Endo-CSV Handpiece, SafeSider HF Instruments in 25 mm length (2 Glidepath Instruments, 2 Shapers, 2 Finishers) 1 box 25/04 GP points, 1 box 30/04 GP points |
| 5100-00 | Endo-CSV Cordless Endodontic Handpiece with LED Illumination and charging base |

SafeSider HF Rotary Endodontic Instruments:

- | | |
|----------|--|
| 5121-123 | SafeSider HF Instruments in 21 mm length (2 Glidepath Instruments, 2 Shapers, 2 Finishers) |
| 5125-123 | SafeSider HF Instruments in 25 mm length (2 Glidepath Instruments, 2 Shapers, 2 Finishers) |
| 5121-01 | SafeSider HF Instruments: 6- 21 mm Glidepath Instruments (20/02) |
| 5125-01 | SafeSider HF Instruments: 6- 25 mm Glidepath Instruments (20/02) |
| 5121-02 | SafeSider HF Instruments: 6- 21 mm Shaper Instruments (30/02) |
| 5125-02 | SafeSider HF Instruments: 6- 25 mm Shaper Instruments (30/02) |
| 5121-03 | SafeSider HF Instruments: 6- 21 mm Finisher Instruments (30/04) |
| 5125-03 | SafeSider HF Instruments: 6- 25 mm Finisher Instruments (30/04) |
| 5104-G25 | SafeSider HF Gutta Percha Points. 60 points. Size 25/04 |
| 5104-G30 | SafeSider HF Gutta Percha Points. 60 points. Size 30/04 |

For our complete line of patented dental products please visit our website at: www.edsdental.com